

Kontakt™



Choose adaptability



BIOTECH DENTAL

THE ICONIC IMPLANT FROM BIOTECH DENTAL

For more than a decade, the Kontakt™ implant has become **Biotech Dental's iconic implant** among dentists, thanks to its **technical** and **clinical performance**, but also thanks to its **adaptability**.



✓ Reliable connection - Morse taper connection (10°)

The Kontakt™ implant is designed with a **conical connection**. Thanks to its **large contact surface**, it ensures the **stability of prosthetic abutments**, prevents micromovements and provides a bacterial seal. The mechanical strength of the implant and abutment is improved.

✓ Patented STSystem® indexing

This indexing system is the **exclusive** property of Biotech Dental. It provides easy repositioning and maximum flexibility thanks to the **6 possible positions** for the prosthetic parts.

✓ Efficient design

The Kontakt™ is a **cylindro-conical** implant with a tapered thread base with compressive targeting to ensure **primary stability** as well as **osseointegration**.

- **Narrowed, chamfered and micro-structured implant neck**
The narrowing neck reduces bone compression of the cortical bone.
The structured chamfer promotes the retention of blood coagulum which is necessary for bone reconstruction.
- **Surface Increaser**
The secondary thread increases the developed surface and promotes the distribution of forces on the bone.
- **Constant Leaf**
Cutting spirals all along the implant improve stability, reducing the effort of insertion into the bone.

✓ Platform switching

Platform switching, in synergy with a stable and bacteria-proof **morse taper** connection, is an important factor concerning the tissue stability. It helps to prevent peri-implantitis.

✓ Atraumatic spherical apex

Thanks to its design, the apex reduces the anatomical risk during surgery and allows a less invasive osteotomy.

2

MADE IN FRANCE WITH A LIFETIME GUARANTEE

Our implants are manufactured at our **Biotech Dental Manufacturing** factory in the Arve Valley, in Scionzier, France.

Our production unit has been maintaining its unique expertise in micro-mechanical manufacturing for 120 years, and our implants are produced with the **rigour, precision and passion inherited from our clockmaking traditions.**



Scionzier



FRENCH ORIGIN GUARANTEED

Our **certification** to the **Origine France Garantie*** label ensures the origin of our medical devices



Lifetime GUARANTEE

Kontakt™ implants are **guaranteed for life** when combined with the use of **original Biotech Dental prosthetic parts** or **Biotech Dental Digital CAD-CAM solutions**



3

BENEFITS

The Kontakt™ is a cylindro-conical implant with primary stability promoting immediate loading, thanks to its design and placement protocol.



IT HAS THE FOLLOWING BENEFITS:

- ✓ **Tissue stability** thanks to its morse taper connection, platform switching and the concave profile of the transgingival part of the abutments. The STSystem® connection prevents bacterial infiltration, preserving the peri-implant biological space.
- ✓ **High mechanical resistance** thanks to the patented design of the STSystem® connection, which creates a large friction surface between the implant and abutment, reducing micromovements. Abutment retention and stability are ensured by the friction between the abutment and the implant.
- ✓ **Prosthetic simplicity and adaptability:** a single STSystem® connection for all implant diameters* simplifies stock management and prosthetic procedures both at the chairside and in the laboratory.

*Except 3mm diameter

4

IMPLANT SCANNING IN YOUR DIGITAL FLOW

- ✓ **Optimise** your digital workflow, communication with your patients and **accelerate** your cases **planning**.
- ✓ **Simplify** your impressions with an intraoral scanner.


With the scanbody you digitally select the implant position and orientation.

Realise perfectly adapted abutments and prosthetic parts to your patient's morphology.



GUIDED SURGERY FOR EVERYONE

Biotech Dental combines the performance of the Kontakt™ implant with NemoScan technology, its proprietary implant planning software, **to give you easy access to guided surgery**.



PLANNING CENTER

An ergonomic, intuitive platform for ordering all types of guides.

www.biotech-planningcenter.com

The planning centre can be used to produce pilot guides, «fully guided» guides and stackable guides.

5

A **SINGLE** SURGICAL AND GUIDED SURGERY KIT FOR ALL IMPLANTS IN THE KONTACT™ RANGE



THE KONTACT™ SURGICAL AND ATLASURGERY™ 2 GUIDED SURGERY KITS ARE:

- ✓ **Compact** for minimum space requirements and optimal storage.
- ✓ **Convenient**, with quick and easy opening for easy access to instruments.
- ✓ **Understandable** thanks to colour marking for quick instrument identification.
- ✓ **Fully removable** for complete decontamination and autoclavability.

6

PROSTHETIC SOLUTIONS FOR ALL INDICATIONS

A COMPLETE PROSTHETIC RANGE DEVELOPED TO SIMPLIFY AND OPTIMISE YOUR AESTHETIC AND FUNCTIONAL RESULTS

- ✓ Cemented and screw-retained prosthesis
- ✓ Removable prosthesis
- ✓ Non-removable prosthesis
- ✓ Anatomical, customisable and scannable SSA-GF* healing abutments
- ✓ Customised CAD-CAM prosthesis



* Sealing Socket Abutment - Gingival Fit



BIOTECH DENTAL DIGITAL, INDUSTRIAL CAD-CAM CENTRE IN THE ARVE VALLEY, IN SCIONZIER - FRANCE 
(Certified ISO13485 and Origine France Garantie*)



- ✓ Single and multiple implant-supported prostheses in Titanium, Cobalt-Chromium and Zirconia
- ✓ Scan design service to support you and your laboratory in processing your digital work
- ✓ Digital prostheses as part of the Biotech Dental guarantee



*French Origin Guaranteed

KONTACT™ A COMPLETE RANGE OF IMPLANTS TO TREAT ALL YOUR CASES


✓ 5 body diameter:



✓ 9 lengths

- 6 mm
- 8 mm
- 10 mm
- 12 mm
- 14 mm
- 16 mm
- 20 mm*
- 22 mm*
- 24 mm*

* Kontakt™ XL implant (20 mm, 22 mm, 24 mm).



XL*

	References	Designations	Diameters	Lengths
	K30-10	Implants	Ø 3 mm	10 mm
	K30-12			12 mm
	K30-14			14 mm
	K3608	Implants	Ø 3,6 mm	8 mm
	K3610			10 mm
	K3612			12 mm
	K3614			14 mm
	K3616			16 mm
	K4206	Implants	Ø 4.2 mm	6 mm
	K4208			8 mm
	K4210			10 mm
	K4212			12 mm
	K4214			14 mm
	K4216*			16 mm
	K4806	Implants	Ø 4.8 mm	6 mm
	K4808			8 mm
	K4810			10 mm
	K4812			12 mm
	K4814			14 mm
	K5406	Implants	Ø 5.4 mm	6 mm
	K5408			8 mm
	K5410			10 mm
	K5412			12 mm
	K5414			14 mm
	K3620	Implants	Ø 3.6 mm	20 mm
	K3622			22 mm
	K3624			24 mm
	K4220	Implants	Ø 4.2 mm	20 mm
	K4222			22 mm
	K4224			24 mm

8 CLINICAL CASES

1 Extractions and immediate loading: Kontakt™ implant in the aesthetic sector


Dr. Elias Khoury



Initial clinical situation



Coronal section highlighting a bone defect



Alveolar curettage and apical anchoring of the implant, alveolar filling and x-ray control




Customised abutments



Aesthetic result on the day of crown placement



Axial cone beam section highlighting the new bone outlines created




X-ray follow-up after 13 years: stability of the regenerated vestibular volume is confirmed



Tissue stability after 13 years

2 Implant rehabilitation in the aesthetic sector


Dr. Guerino Paolantoni



Initial case: 35 year-old woman
Left upper incisor with vertical fracture



Implant surgery: after removing the complete flap, the bone defect revealed a vestibular dehiscence of 15 mm




Implant placement: use of a healing screw to improve the stability of the blood clot



Guided bone regeneration: use of a resorbable collagen membrane to stabilise the bovine xenograft



The temporary screwed prosthesis was placed six months after the implant surgery, in order to assist the healing of the peri-implant tissues



View of healed peri-implant tissues




Final screw-retained crown and x-ray follow-up after 7 years

3

SSA-GF* workflow: Step-by-step treatment of a mandibular molar replacement


Dr. Gary Finelle




Initial situation: 46 cannot be preserved. Distal root fracture




Atraumatic extraction. Root separation



Intraseptal implant drilling. Final drilling



Placing the implant in the ideal prosthetic position




Fitting the mandibular SSA-GF abutment. The abutment seals almost the entire alveolar emergence



Customising the SSA-GF abutment with a flow composite to seal



Final placement of the SSA-GF abutment after customisation: mechanical barrier between alveolar site and oral cavity, tissue support, blood clot stabilisation



Placement of the monolithic zirconia implant crown

4

Immediate loading and alveolar preservation around a Kontakt™ implant

Dr. Hadi Antoun



36 is in the terminal phase, an immediate loading is scheduled in order to reduce the number of procedures, shorten the overall treatment time and maintain bone and mucosal volume as much as possible




A Kontakt™ implant (diameter 4.2 mm - length 10 mm) is placed after the least traumatic extraction possible and careful cleaning of the alveolus. Filling and placement of a PEEK SSA-GF abutment



At 4 months, the implant is inspected for healing and integration. A digital impression is then taken by the practitioner, directly on the SSA-GF abutment to record the position of the implant and the anatomical emergence profile



A stratified zirconia crown is made, bonded to a Ti-base and screwed directly onto the implant




Retro-alveolar follow-up after one year to check the x-ray integration of the Kontakt™ implant as well as the stability of the stabilised proximal bone at and above the implant neck (Dr. Frédéric Pujol Prosthesis)

5

FP1 rehabilitation: Socket Shield technique in a patient with paradontal disease


Dr. Helder Moura



6

Full rehabilitation

Dr. Bertrand Rousselet




7 | Stackable guide - YOUR3DGUIDE™

Dr. Pierre Keller



Preoperative situation



Planning performed by the Biotech Dental planning centre via NemoScan




Placement of the stackable guide's titanium base with the positioning stage



Placement of Kontakt™ implants with the drilling stage and appropriate implant holders



Placement of conical abutments




Connection of the temporary prosthesis (last stage of the guide) on conical abutments via the sleeves



Finishing the machined PMMA temporary prosthesis




Immediate loading




Panoramic x-ray after placing Kontakt™ implants

8 | YOUR3DCAGE™ bone reconstruction and implant placement

Dr. Pierre Keller




Pre-operative cone-beam assessment




Bone volume design




Design of the Your3Dcage™ titanium membrane




Careful preparation of the bone defect



Filling the titanium membrane with autogenous bone



Fixing the Your3Dcage™ titanium membrane with 2 screwpins




Removal of Your3Dcage™ membrane after 4 months and visualisation of bone gain




Placement of the pilot guide




Placing the first Kontakt™ implant



Placing 2 high cover screws on Kontakt™ implants



Guided bone regeneration using Nea Cova™ membrane fixed by 2 screwpins



Post-operative panoramic x-ray



BIOTECH DENTAL



305, Allées de Craponne
13300 Salon-de-Provence - France
Tel.: +33 (0)4 90 44 60 60
Fax: +33 (0)4 90 44 60 61

CustomerCareService@biotech-dental.com
www.biotech-dental.com

www.biotechgalaxy.com



AFNOR Cer 73017

Biotech Dental - 305, Allées de Craponne - 13300 Salon-de-Provence, France. A simplified company limited by shares (S.A.S.) with a capital of 24,866,417 € - Corporate and Trade Register of Salon de Provence: 795 001 304 - SIRET: 795 001 304 00018 - VAT no.: FR 31 79 500 13 04.

Kontakt™: manufactured by Biotech Dental. Class I, IIa and IIb medical devices. CE0459. NemoScan: manufactured by Software Nemotec, S.L. Class Im MD. CE2797. All of these medical devices must be used by qualified and trained health professionals. All brands are under the responsibility of their respective manufacturers. Read the instructions carefully. No refund from French social security.

Biotech Dental Digital - 305 Allées de Craponne - 13300 Salon de Provence - FRANCE

Tel.: +33 (0)4 65 84 46 56 - info@biotech-dental-digital.com - www.biotech-dental-digital.com

A simplified company limited by shares (S.A.S.) with a capital of 10,000 € - Corporate and Trade Register of Salon de Provence: 797 755 188 - SIRET: 797 755 188 00043 - VAT no.: FR 96 797 755 188

Non contractual visuals. Do not litter in public areas. Printed by VALLIERE - 163, Avenue du Luxembourg - ZAC des Molières - 12140 Miramas - France.