

4 STEPS

to regain your smile

- 1.** To start with, your practitioner will be **taking impressions** of your teeth.
- 2.** Then, he/she will propose a **treatment plan** adapted to your case to **create your aligners**. You will also have access to your smile before/after video.
- 3.** Your treatment will last between **6 and 24 monts**. Aligners should be worn **22 h/per day and during 2 to 3 weeks**. Be careful, you must remove your aligners only when eating and when brushing your teeth.
- 4.** Committing with us means **obtaining and keeping the smile you want to have** with complete discretion. You will need to maintain the treatment effects with **retainers**.

*This video is a treatment simulation. Compliance with the instructions for using the device and compliance with the treatment plan created by the healthcare professional are essential in order to obtain the expected result. Additional phases may be required.

smilers®

Le sourire chic à la française

Biotech Dental Smilers® is a family mid-sized company with strong values, committed to provide a French sourire chic to everyone.

PRACTITIONER'S STAMP



Biotech Dental Smilers
305, Allées de Craponne
13300 Salon-de-Provence
FRANCE

+33 (0)4 13 22 82 40

info@smilers.com



Aligners manufactured by Biotech Dental Smilers. Class IIa custom-made medical devices for dental alignment. Carefully read the instructions on the leaflets. Seek advice from your practitioner. Images are for representation purpose only. Biotech Dental Smilers - Biotech Dental Smilers - 305, Allées de Craponne - 13300 Salon-de-Provence - France. A French company under the french law with a capital of 125 000 € - Commercial register Salon-de-Provence: 801 674 375 - SIRET: 801 674 375 00010 - VAT number: FR 93 80 167 437 5. Do not throw in public areas. Imprimerie VALLIERE - 163, Avenue du Luxembourg - ZAC des Molières - Miramas - France.

SMI-DE-EN_Rév.00_01/2022

smilers®



Regain your
power to smile

Le sourire chic à la française

WHY START A SMILERS® CLEAR ORTHODONTIC TREATMENT ?

Your aligners will help you **realign your teeth** to get the smile you deserve.

These removable and clear aligners will adjust to your teeth **discreetly. They will be custom-made.**

MANY TREATABLE CASES



In some cases, your practitioner may need to:

- > Glue small resin dots called **attachments** to **ease some movements.**
- > Perform a **stripping***: file a few tenths of a millimeter of enamel to allow some movements. **This practice won't hurt you.**

*interproximal reduction.



ECO- RESPONSABILITY

Our aligners are bisphenol-free, their packaging is **recyclable** and our brochures are printed **with vegetable-based ink.**

MADE IN FRANCE QUALITY



✓ **Aligners inspired by french Haute Couture for your comfort:**

- Produced in our **high technology factory** in Salon-de-Provence
- **Precise** cut along the dental neck
- Hand-polished in our Smilers® workshop

✓ **A smooth impression which prevents deposits and preserves the clearness day after day**

Our factory is **certified ISO 13485**, proving the qualitative production of our clear aligners



Smilers® aligners compared to competition aligners, after one week in coffee.