

SSA-GF*

custom-made healing



The **only anatomical, customisable** and **scannable** healing abutments solution.

ANATOMICAL¹

CUSTOMISABLE



SCANNABLE²

COMPATIBLE

with Kontakt™, Nobel RP, Straumann RC, Straumann RN, Straumann WN.

*Sealing Socket Abutment - Gingival Fit.

¹ Solution for immediate post-extraction implant placement: no incision, no stitches, no membrane.

² Prosthetic simplification thanks to scannable caps. Exempt from references compatible with Straumann and Nobel implants.



Kontakt™ «direct implant» parts and Nobel/Straumann compatibility are not available in all countries. Please contact your local representative for further information.



BIOTECH DENTAL


SOLUTION « ENDO-OSSEOUS IMPLANTS » KONTACT™ RANGE

SSA-GF* DIRECT IMPLANT



- ✓ Biocompatible material (PEEK)
- ✓ Anatomic
- ✓ Personalised
- ✓ Scannable
- ✓ Tight seal morse-taper locking system (ST System®)
- ✓ Compatible with **Kontakt™ Kontakt™ N, Kontakt™ S, Kontakt™ S+**
- ✓ Recommended in molar and pre-molar areas
- ✓ Recommended for immediate post-extraction implant placement

SSA-GF* ON OMNIPOST ABUTMENT



- ✓ One Abutment - One Time
- ✓ Biocompatible material (PEEK)
- ✓ Anatomic
- ✓ Personalised
- ✓ Scannable
- ✓ Compatible with **Kontakt™ Kontakt™ N, Kontakt™ S, Kontakt™ S+**
- ✓ Recommended in molar and pre-molar areas
- ✓ Recommended for immediate post-extraction implant placement

SSA-GF* Direct implant Kontakt™

	References	Designations	Shapes	Heights
	BK70-4H2	Healing screw BK 4H2	Universal	2 mm
	BK70-5H2	Healing screw BK 5H2	Mandibular	2 mm
	BK70-65H2	Healing screw BK 6,5H2	Maxillary	2 mm
	BK70-4H4	Healing screw BK 4H4	Universal	4 mm
	BK70-5H4	Healing screw BK 5H4	Mandibular	4 mm
	BK70-65H4	Healing screw BK 6,5H4	Maxillary	4 mm

SSA-GF* anatomic healing caps

	References	Descriptions
	KPOPSSA	SSA-GF cap for OmniPost abutment
	KPOPSSAMD	SSA-GF mandibular molar cap for OmniPost abutment
	KPOPSSAMX	SSA-GF maxillary molar cap for OmniPost abutment

Straight indexed OmniPost abutments

	References	Designations	Angulations	Heights
	KPOP150	Straight indexed OmniPost abutments + KVP screw	0°	1.5 mm
	KPOP200			2 mm
	KPOP300			3 mm
	KPOP400			4 mm
	KPOP500			5 mm

SOLUTION « KONTACT™ PERIO LEVEL » SSA-GF* DIRECT IMPLANT




- ✓ Biocompatible material (PEEK)
- ✓ « Tissue-Level » implant solution
- ✓ Anatomic
- ✓ Personalised
- ✓ Several gum heights available
- ✓ Scannable
- ✓ Recommended in molar and pre-molar areas
- ✓ Recommended for immediate post-extraction implant placement

SSA-GF* abutments for Kontakt™ Perio Level implant

	References	Designations	Neck diameters	Heights
Universal	KLPSSA-S	SSA-GF abutment + KPLVP screw	S	2 mm
	KLPSSA-M	SSA-GF abutment + KPLVP screw	M	2 mm
Mandibular	KLPSSAMD-L	SSA-GF mandibular molar abutment + KPLVP screw	L	1 mm
	KLPSSAMD2-L	SSA-GF mandibular molar abutment neck L Ht. 2 mm + KPLVP Screw		2 mm
	KLPSSAMD-XL	SSA-GF mandibular molar abutment + KPLVP screw	XL	1 mm
	KLPSSAMD2-XL	SSA-GF mandibular molar abutment neck XL Ht. 2 mm + KPLVP screw		2 mm
Maxillary	KLPSSAMX-L	SSA-GF maxillary molar abutment + KPLVP screw	L	1 mm
	KLPSSAMX2-L	SSA-GF maxillary molar abutment neck L Ht. 2 mm + KPLVP screw		2 mm
	KLPSSAMX-XL	SSA-GF maxillary molar abutment + KPLVP screw	XL	1 mm
	KLPSSAMX2-XL	SSA-GF maxillary molar abutment neck XL Ht. 2 mm + KPLVP screw		2 mm

*Sealing Socket Abutment - Gingival Fit.

NOBEL & STRAUMANN COMPATIBILITY



- ✓ Biocompatible material (PEEK)
- ✓ Anatomic
- ✓ Personalised
- ✓ Recommended in molar and pre-molar areas
- ✓ Recommended for immediate post-extraction implant placement

NOBEL COMPATIBILITY

NOBEL RP

MANDIBULAR MOLAR



Ht. 1 mm



Ht. 3 mm



MAXILLARY MOLAR



SSA-GF* Compatibility with Nobel Active RP

	References	Designations	Implants type	Shapes	Heights
	NB70-RP-5H2	Healing screw NB RP 5H2	Nobel RP	Mandibular	1 mm
	NB70-RP-5H5	Healing screw NB RP 5H5	Nobel RP	Maxillary	1 mm
	NB70-RP-5H3	Healing screw NB RP 5H3	Nobel RP	Mandibular	3 mm
	NB70-RP-5H6	Healing screw NB RP 5H6	Nobel RP	Maxillary	3 mm

STRAUMANN COMPATIBILITY

BONE LEVEL

STRAUMANN RC

MANDIBULAR MOLAR



Ht. 1 mm



Ht. 3 mm



MAXILLARY MOLAR



TISSUE LEVEL

STRAUMANN RN

MANDIBULAR MOLAR



Ht. 2 mm



MAXILLARY MOLAR




STRAUMANN WN


MANDIBULAR MOLAR



Ht. 1 mm



MAXILLARY MOLAR



SSA-GF* Compatibility with Straumann

	References	Designations	Implants type	Shapes	Heights
	SB70-RC-6H2	Healing screw SB RC 6H2	Straumann RC	Mandibular	1 mm
	SB70-RC-6H4	Healing screw SB RC 6H4	Straumann RC	Maxillary	1 mm
	SB70-RC-6H3	Healing screw SB RC 6H3	Straumann RC	Mandibular	3 mm
	SB70-RC-6H5	Healing screw SB RC 6H5	Straumann RC	Maxillary	3 mm
	ST70-RPH2	Healing screw ST RP H2	Straumann RN	Mandibular	2 mm
	ST70-RPH3	Healing screw ST RP H3	Straumann RN	Maxillary	2 mm
	ST70-WPH2	Healing screw ST WP H2	Straumann WN	Mandibular	1 mm
	ST70-WPH3	Healing screw ST WP H3	Straumann WN	Maxillary	1 mm