


IRIS BY STARCK®



Innovation. Design. Technology.

Design & technology



CEO of Biotech Dental Group

- Philippe Veran

The best way to predict the future is to shape it"



Creator

- Philippe Starck


The spirit is my favourite means of transport, because it is unlimited and the battery is never empty"

“Invented, designed,
and manufactured
in France

Philippe Starck

“French engineering at its finest”


Philippe Starck



Optical design

Design & Ergonomics

Power meets elegance



Biomimetic curve:
handpiece that is an extension of the hand

Fin: better grip and easier handling for greater precision in movement

Remote actions:
minimise interruptions for a smoother workflow

LED indicator light:
indicates whether the camera is in standby mode or capturing images

Fully transparent case:
reveals the finesse of high-tech components

Rounded outlines:
no sharp angles for a softer feeling in the mouth and a better patient experience

“organic design”

Philippe Starck

RIS by STARCK

Controlled environmental impact

- Hosting management on internally controlled servers
- No disposable tips: sterilisation by immersion (IP67 sealing) and UVC to reduce waste during use (no additional costs)
- Considering environmental issues
- Designed without batteries: continuous power source without creating waste that is difficult to recycle

“In perfect harmony with the body”

Philippe Starck



Artificial intelligence & Scalable technology

Artificial intelligence



- Automatic detection of shapes and textures (especially scanbodies)
- High scanning speed
- Ultimate detail identification
- Ability to map an environment in real time (the camera knows its exact position in real time in the mouth)
- High precision of 3D models
- Scalable device: Immediate access to technological innovations and AI-enhanced features
- OLI: provides Level I technical support

“A visionary innovation”

Philippe Starck

Security & Interoperability

- **Secure data sharing environment**
- **Customised team: adding secure accounts for colleagues and assistants**
- **Scalable, highly secure storage (Scan4All) in France**
- **ISO 13485, ISO 27001, GDPR, HDS certifications**


“serving both
patient
and dentist
with care

Philippe Starck




“Where technology,
ergonomics &
humanity meet

Philippe Starck



An ergonomic ecosystem



IRIS, A true revelation



World-unique patented technologies



Made in France



Designed by Philippe Starck



International communication platform: IRIS Connect



Specifications

CATEGORY	IRIS BY STARCK
CONFIGURATION	POD (Kart version, coming soon)
HANDPIECE WEIGHT	120 g
HANDPIECE SHAPE TRANSPARENT FRONT PART DIMENSION (L x l x H)	Organic shape designed by Starck 6 * 22 * 16 mm
CONNECTION PORT	USB-C
WIRELESS	NO
TIP	NO
DISINFECTION	SOAKING (UV-C available soon) IP67 SEALING
RESOLUTION	1200 * 1200 px (HD sensors)
depth of field	20 mm
CALIBRATION	NO
SCAN PRINCIPLE	PASSIVE IMAGING
DATA OUTPUT	STL, PLY
SOFTWARE	EASY & INTUITIVE
SECURE DATA EXCHANGE	YES (ISO 27 001, HDS) Free and open to all your contacts
MULTI USER	YES, with right management
ARTIFICIAL INTELLIGENCE	YES
GUARANTEE	1 YEAR (warranty extension available)

IRIS BY STARCK®

IRIS: manufacturer Scan4all SAS. Class I medical device dedicated to digital dentistry. CE.
Scan4all - 137 rue Garance - 31670 Labège, France. A simplified company limited by shares (S.A.S.) with a capital of 191,590 € -
SIRET: 50748486300059 - VAT No.: FR25507484863 - Corporate and Trade Register of Toulouse. Non contractual visuals. Do not litter in public areas.
Printed by VALLIERE - 163, Avenue du Luxembourg - ZAC des Molières - 13140 Miramas - France.