


SIMPLIFYING YOUR CLINICAL PROTOCOLS

BioscanHealer eliminates the need to remove the healing abutment, simplifying the procedure and saving valuable time.



1 Place



2 Scan



3 Restore



BIOTECH DENTAL
GROUP

Biotech Dental
305, Allées de Craponne
13300 Salon-de-Provence - France


Tel.: +33 (0)4 90 44 60 60 | Fax: +33 (0)4 90 44 60 61

CustomerService@biotech-dental.com
www.biotech-dental.com
www.galaxybiotech.com

BIOTECH DENTAL
IMPLANTS

X

BIOSCAN
HEALER



BIOSCAN
HEALER

SYNCHRONISING HEALING
WITH THE POWER OF DIGITAL



BIOTECH DENTAL
GROUP

CREATING A NATURAL TISSUE EMERGENCE

The protocol has been simplified and the patient process considerably streamlined for a more comfortable and efficient experience.


INTEGRATED DIGITISATION

BioscanHealer has been designed to provide optimal accuracy.


6 available heights*

Easily identified thanks to a dot system on the scannable side.

*According to diameters



Points should be positioned on the palatal side.



4 available diameters

Easily identified thanks to the different variations of scannable flat surfaces.

FULL RANGE



Available libraries
info@biotech-dental-digital.com
 +33 (0)4 13 22 16 29

